



5

IHE <Domain Name> Technical Framework Supplement

10

<Profile Name
(Profile Acronym) >

<For FHIR based profiles, indicate the FHIR release number & the FMM levels of the contents; otherwise, delete the two following lines.>

15

HL7[®] FHIR[®] STU **x**

Using Resources at FMM Level **n-n**

Revision x.x – Draft in Preparation for Public Comment (*or* Trial Implementation)

20

<The IHE Documentation Specialist will change the title to just “Draft for Public Comment” or “Trial Implementation” upon publication. Leave “as is” until then.>

Date: <Month xx, 20xx>

Author: <Domain> Technical Committee

Email: <domain_name@ihe.net>

25

Please verify you have the most recent version of this document. See [here](#) for Trial Implementation and Final Text versions and [here](#) for Public Comment versions.

<Instructions to authors are encapsulated in angled brackets as “< ... >” and denoted with italicized text. These instructions should be deleted entirely prior to publication.>

- 30 *<Use of capitalization: Please follow standard English grammar rules-only proper nouns and names are upper case. For example, “Modality Actor” is upper case, but “an actor which fulfills the role of a modality” is lower case. Do not use upper case to emphasize a word/topic. Examples:*

- 35 *<Note: Before creating a draft supplement, please review the editing conventions, which include information such as section, table and diagram numbering and how to use Microsoft Word tools, at http://wiki.ihe.net/index.php?title=Writing_Technical_Frameworks_and_Supplements. This guidance is especially useful for first time authors.>*

- 40 *<This supplement template is intended for developing new profiles or making significant changes to profiles, such as adding formal options. Simple changes to existing supplements or profiles should be made using the Change Proposal (CP) process. See the Technical Framework Development section at http://wiki.ihe.net/index.php?title=Process#Technical_Framework_Development for more guidance on supplements vs. CPs.>*

- 45 *<All of the sections in this document are required. Sections may not be deleted. The outline numbering is intended to be consistent across profiles and across domains, so do not adjust the outline numbering. If there is no relevant content for a section, simply state “Section not applicable”, but leave the numbering intact. Sub-sections may be added for clarity.>*

<This supplement template includes templates for Volumes 1 (Profiles), 2 (Transactions), 3 (Content Modules), and 4 (National Extensions).>

- 50 *<Volumes 1, 2, and/or 3 are developed together for Public Comment and Trial Implementation submission. Volume 4, National Extensions, is typically developed at a later point in time, usually at Trial Implementation or later. Templates for all four volumes are included in this document for the sake of completeness. If you are beginning a new profile, you are strongly discouraged from using National Extensions and should instead focus on optional data sets or*
- 55 *other alternatives. For more information, see http://wiki.ihe.net/index.php?title=National_Extensions_Process.>*

Foreword

This is a supplement to the IHE <Domain Name> Technical Framework <VX.X>. Each supplement undergoes a process of public comment and trial implementation before being incorporated into the volumes of the Technical Frameworks.

<For Public Comment:> This supplement is published on <Month XX, 201x> for Public Comment. Comments are invited and can be submitted at http://www.ihe.net/Public_Comment/#domainname. In order to be considered in development of the Trial Implementation version of the supplement, comments must be received by <Month XX, 201X>.

<For Trial Implementation:> This supplement is published on <Month XX, 201X> for Trial Implementation and may be available for testing at subsequent IHE Connectathons. The supplement may be amended based on the results of testing. Following successful testing it will be incorporated into the <Domain Name> Technical Framework. Comments are invited and can be submitted at http://www.ihe.net/Public_Comment/#domainname.

This supplement describes changes to the existing technical framework documents.

“Boxed” instructions like the sample below indicate to the Volume Editor how to integrate the relevant section(s) into the relevant Technical Framework volume.

Amend section X.X by the following:

Where the amendment adds text, make the added text **bold underline**. Where the amendment removes text, make the removed text **~~bold strikethrough~~**. When entire new sections are added, introduce with editor’s instructions to “add new text” or similar, which for readability are not bolded or underlined.

General information about IHE can be found at www.ihe.net.

Information about the IHE <Domain Name> domain can be found at ihe.net/IHE_Domains.

Information about the organization of IHE Technical Frameworks and Supplements and the process used to create them can be found at http://ihe.net/IHE_Process and <http://ihe.net/Profiles>.

The current version of the IHE <Domain name> Technical Framework can be found at http://ihe.net/Technical_Frameworks.

<Comments may be submitted on IHE Technical Framework templates any time at http://ihe.net/Templates_Public_Comments. Please enter comments/issues as soon as they are found. Do not wait until a future review cycle is announced.>

CONTENTS

	Introduction to this Supplement.....	7
95	Volume 1 – Profiles	8
	X.4.2 Use Cases	8
	X.4.2.1 Use Case #1: SDPi-R Standalone OR Dashboard Use Case - SORD	8
	X.4.2.1.1 SDPi SORD Use Case Description.....	8
	X.4.2.1.2 SDPi SORD Process Flow	8
100	X.4.2.2 Use Case #2: SDPi-R Standalone OR Cockpit Use Case - SORC	10
	X.4.2.2.1 SDPi SORC Use Case Description.....	10
	X.4.2.2.2 SDPi SORC Process Flow	11
	X.4.2.3 Use Case #3: SDPi-R/xC ICU Isolation PoC Use Case - IIPoC	14
	X.4.2.3.1 SDPi IIPoC Use Case Description.....	14
105	X.4.2.3.2 SDPi IIPoC Process Flow	14
	X.4.2.4 Use Case #4: SDPi-A Silent PoC Alert Distribution Use Case - SPoC	17
	X.4.2.4.1 SDPi SPoC Use Case Description	17
	X.4.2.4.2 SDPi SPoC Process Flow	17
	X.4.2.5 Use Case #5: SDPi-A Silent ICU Alert Distribution - SICU	20
110	X.4.2.5.1 SDPi SICU Use Case Description	20
	X.4.2.5.2 SDPi SICU Process Flow	20
	X.4.2.6 Use Case #6: SDPi-A/xC Remote Alert Management - RAM.....	23
	X.4.2.6.1 SDPi RAM Use Case Description	23
	X.4.2.6.2 SDPi RAM Process Flow	23
115	X.4.2.7 Use Case #7: SDPi-A Smart Alerting System - SAS	25
	X.4.2.7.1 SDPi SAS Use Case Description	25
	X.4.2.7.2 SDPi SAS Process Flow	25
	Appendices to Volume 1	30
	Appendix C – Use Cases - Detail	31
120	C.1 Overview of the Concepts for DIS and DAS	31
	C.2 Feature: Use Case 1 - SDPi-R Standalone OR Dashboard - SORD.....	34
	C.2.1 Narrative:	34
	C.2.2 Background: Technical Pre-Conditions.....	34
	C.2.3 Scenario: SORD 1.1 - OR Devices are Accessible to the Dashboard	34
125	C.2.4 Scenario: SORD 1.2 - OR Devices are inaccessible to the Dashboard	34
	C.2.5 Scenario: SORD 1.3 - One or more OR Devices are inaccessible to the Dashboard	34
	C.3 Feature: Use Case 2 - SDPi-R/A/xC Standalone OR Cockpit - SORC.....	34
	C.3.1 Narrative:	35
	C.3.2 Background: Technical Pre-Condition	35
130	C.3.3 Scenario: SORC 2.1 - OR Devices are accessible to the Cockpit	35
	C.3.4 Scenario: SORC 2.2 - OR Devices are inaccessible to the Cockpit	36
	C.3.5 Scenario: SORC 2.3 - Some OR Devices are inaccessible to the Cockpit	36
	C.4 Feature: Use Case 3 - SDPi-R/xC ICU Isolation PoC Use Case - IIPoC.....	36
	Narrative:.....	36

135	Background: Technical Pre-Condition.....	36
	C.5 - Feature: Use Case 4 - SDPI-A - Silent PoC Alert distribution - SPoC.....	36
	C.5.1 Narrative	37
	C.5.2 Background: Technical Pre-Condition	37
140	C.5.3 Scenario: SPoC 4.1 - Medical device detects an alert situation and at least one distributed alarm system (xDASxx) is accessible	38
	C.5.4 Scenario: SPoC 4.2 - Medical device detects an alert situation and all distributed alarm systems (xDASxx) are inaccessible or become inaccessible.....	38
	C.5.5 Scenario: SPoC 4.3 - Alert situation has been resolved and at least one distributed alarm system (xDASxx) is accessible.....	38
145	C.5.6 Scenario: SPoC 4.4 - Medical device detects an alert situation, initially xDASxx is accessible but fails	39
	C.6 Feature: Use Case 5 - SDPi-A: Silent ICU Alert Distribution - SICU	39
	C.6.1 Narrative:	39
	C.6.2 Background: Technical Pre-Condition	40
150	C.6.3 Scenario: SICU 5.1 - Medical device detects an alert situation and the distributed alarm system CDAScg is accessible.....	40
	C.6.4 Scenario: SICU 5.2 - Medical device detects an alert situation and the distributed alarm system CDAScg is inaccessible.....	40
	C.6.5 Scenario: SICU 5.3 - Medical device detects an alert situation and all distributed alarm systems (xDASxx) are or become inaccessible	40
155	C.6.6 Scenario: SICU 5.4 - Alert situation has been resolved and at least one distributed alarm system (xDASxx) is accessible.....	41
	C.6.7 Scenario: SICU 5.5 - Medical device detects an alert situation, initially DAS is accessible but fails	41
160	C.7 Feature: Use Case 6 - SDPi-A/xC: Remote Alert Management - RAM.....	41
	C.7.1 Narrative	41
	C.7.2 Background: Technical Pre-Condition	42
	C.7.3 Scenario: RAM 6.1 - Caregiver adjusts alarm limit at their Mobile Device	42
165	C.7.4 Scenario: RAM 6.2 - Caregiver attempts to adjust alarm limit at their CDAScg mobile device but the respective CDAScg is not accessible	43
	C.8 Feature: Use Case 7 - SDPi-A Smart Alerting System - SAS	43
	C.8.1 Narrative:	43
	C.8.2 Background: Technical Pre-Conditions.....	44
170	C.8.3 Scenario Outline: SAS 7.1 - Local device generates alerts, Smart Alerting System is accessible, and xDASmp and/or xDAScg is accessible	44
	C.8.3 Scenario Outline: SAS 7.2 - Local device generates alerts, Smart Alerting System is inaccessible, and xDASmp and/or xDAScg is accessible	45
	C.8.4 Scenario: SAS 7.3 - Local device generates alerts, and xDASmp and/or xDAScg is inaccessible	45
175	C.8.5 Scenario: SAS 7.4 -Smart Alerting generates an alert and an xDASxx is accessible	46
	C.8.6 Scenario: SAS 7.5 - Smart Alerting generates an alert and an xDASxx is inaccessible	

180	C.8.6 Scenario: SAS 7.6 - Smart Alerting System is inaccessible from devices and an xDASxx is accessible	46
	C.8.7 Scenario: SAS 7.7 - Smart Alerting System is inaccessible from devices and an xDASxx is inaccessible	46
	C.8.8 SAS Scenario 7.8: Devices, CDASxx and Smart Alerting System are accessible, and SAS generates an alert	47

Introduction to this Supplement

185

Volume 1 – Profiles

X.4.2 Use Cases

For detailed information concerning the SDPi Use Cases, Narratives and Scenarios please refer to Volume 1 - Appendix C.

190 X.4.2.1 Use Case #1: SDPi-R Standalone OR Dashboard Use Case - SORD

Aggregated display of information relevant to a surgeon in an Operating Room from multiple data sources on a “Dashboard”.

Additional information on this Use Case can be found in Volume 1 – Appendix C.2.

X.4.2.1.1 SDPi SORD Use Case Description

195 This Use Case describes a general purpose “Dashboard” information display which aggregates information from different sources. While the Dashboard may display visual alert information it does not signal alerts or have any remote-control capabilities.

X.4.2.1.2 SDPi SORD Process Flow

200 *<Diagram and describe the process flow(s) covered by this profile in order to satisfy the use cases. Demonstrate how the profile transactions are combined/sequenced. To provide context and demonstrate how the profile interacts with other profiles, feel free to include transactions and events that are “external” to this profile (using appropriate notation.)*

205 *The set of process flows will typically be exemplary, not exhaustive (i.e., it will address all the use cases, but will not show all possible combinations of actors, or all possible sequencing of transactions).*

If there are detailed behavioral rules that apply to a specific process flow or multiple process flows, an appendix may be added as needed.>

<The roles at the top of the swimlane diagram should correspond to actor names, include the profile acronym:actor name if referencing an actor from a different profile.>

210 *<Modify the following “Swimlane Diagram”.>*




Figure X.4.2.1-1: Basic Process Flow in SDPI Profile – SORD Use Case

215 *<If process flow “swimlane” diagrams require additional explanation to clarify conditional flows, or flow variations need to be described where alternate systems may be playing different actor roles, document those conditional flows here.>*

<Delete the material below if this is a workflow or transport profile. Delete the material above if this profile is a content module only profile.>

220

Pre-conditions:

- Bedside devices are SDC compliant
- At least one OR Dashboard (DISp) (see Annex C)

225 **Main Flow:**

- Devices are all turned on.
 - Devices advertise their presence
- User sets the Dashboard device into configuration mode.
 - Dashboard device builds list of available devices in the room.
- 230 ○ User configures which devices he/she wants to view on the Dashboard.
- User sets the Dashboard device into operational mode.
 - Dashboard connects to each configured device and collects the data from each device and displays it on the Dashboard display.

235 **Post-conditions:**

The surgeon has an overview of all relevant data regarding the patient and device status.

X.4.2.2 Use Case #2: SDPi-R Standalone OR Cockpit Use Case - SORC

240 Aggregated display of information relevant to a surgeon in an Operating Room from multiple data sources on a Surgical “Cockpit”.

Additional information on this Use Case can be found in Volume 1 – Appendix C.3.

X.4.2.2.1 SDPi SORC Use Case Description

245 This Use Case discusses a general purpose Surgical “Cockpit” information and control display which aggregates information from different SDC compliant devices. The Cockpit can signal audio alerts and may remotely control the SDC compliant devices.

X.4.2.2.2 SDPi SORC Process Flow

250 *<Diagram and describe the process flow(s) covered by this profile in order to satisfy the use cases. Demonstrate how the profile transactions are combined/sequenced. To provide context and demonstrate how the profile interacts with other profiles, feel free to include transactions and events that are “external” to this profile (using appropriate notation.)*

The set of process flows will typically be exemplary, not exhaustive (i.e., it will address all the use cases, but will not show all possible combinations of actors, or all possible sequencing of transactions).

255 *If there are detailed behavioral rules that apply to a specific process flow or multiple process flows, an appendix may be added as needed.>*

<The roles at the top of the swimlane diagram should correspond to actor names, include the profile acronym:actor name if referencing an actor from a different profile.>

<Modify the following “Swimlane Diagram”.>

260




Figure X.4.2.2-1: Basic Process Flow in SDPi Profile – SORC Use Case

265 *<If process flow “swimlane” diagrams require additional explanation to clarify conditional flows, or flow variations need to be described where alternate systems may be playing different actor roles, document those conditional flows here.>*

<Delete the material below if this is a workflow or transport profile. Delete the material above if this profile is a content module only profile.>

270 **Pre-conditions:**

- Bedside devices are SDC compliant
- At least one Anesthesia and/or Surgical Cockpit (xDASsp) (see Annex C)

Main Flow:

- 275
- Devices are all turned on.
 - Devices advertise their presence
 - User sets the Cockpit device into configuration mode.
 - Cockpit device builds list of available devices in the room.
 - User configures which devices he/she wants to view on the Cockpit.
- 280
- Cockpit advises viewed devices that it will be alerting on their behalf.
 - User sets the Cockpit device into operational mode.
 - Cockpit connects to each configured device and collects the data from each device and displays it on the Cockpit display.
 - Device goes into alarm
- 285
- If a device viewed by the Cockpit goes into alarm, the Cockpit will signal the alarm.
 - The Cockpit will acknowledge receipt and processing of the alarm back to the alarming device.
 - The alarming device may or may not signal the alert based on its configuration.
- 290
- User decides to Accept the alarm
 - If the user decides to accept the alarm, the Cockpit will send an Accept message to the alarming device and stop signaling the alarm locally.

Post-conditions:

295 The surgeon or anesthesiologist has a broad overview and a single point of control for devices in the room.

X.4.2.3 Use Case #3: SDPi-R/xC ICU Isolation PoC Use Case - IIPoC

300 A remotely located display provides access to data from patient connected devices and remote control of those devices.

Additional information on this Use Case can be found in Volume 1 – Appendix C.4.

X.4.2.3.1 SDPi IIPoC Use Case Description

305 This Use Case covers the general case of a remote Cockpit which allows access to patient data and remote control of a patient’s devices in an isolation situation.

X.4.2.3.2 SDPi IIPoC Process Flow

310 *<Diagram and describe the process flow(s) covered by this profile in order to satisfy the use cases. Demonstrate how the profile transactions are combined/sequenced. To provide context and demonstrate how the profile interacts with other profiles, feel free to include transactions and events that are “external” to this profile (using appropriate notation.)>*

The set of process flows will typically be exemplary, not exhaustive (i.e., it will address all the use cases, but will not show all possible combinations of actors, or all possible sequencing of transactions).

315 *If there are detailed behavioral rules that apply to a specific process flow or multiple process flows, an appendix may be added as needed.>*

<The roles at the top of the swimlane diagram should correspond to actor names, include the profile acronym:actor name if referencing an actor from a different profile.>

<Modify the following “Swimlane Diagram”.>




Figure X.4.2.3-1: Basic Process Flow in SDPi Profile – IIPoC Use Case

325 *<If process flow “swimlane” diagrams require additional explanation to clarify conditional flows, or flow variations need to be described where alternate systems may be playing different actor roles, document those conditional flows here.>*

<Delete the material below if this is a workflow or transport profile. Delete the material above if this profile is a content module only profile.>

Pre-conditions:

- 330
- Bedside devices are SDC compliant
 - At least one xDASp (see Annex C) that is SDC compliant
 - Devices in room allow remote-control from the sDASp.

Main Flow:

- 335
- Devices are all turned on.
 - Devices advertise their presence
 - User sets the Cockpit device into configuration mode.
 - Cockpit device builds list of available devices in the room.
 - User configures which devices he/she wants to view on the Cockpit.
 - 340 ○ Cockpit advises viewed devices that it will be alerting on their behalf.
 - User sets the Cockpit device into operational mode.
 - Cockpit connects to each configured device and collects the data from each device and displays it on the Cockpit display.
 - User decides to adjust settings on a (remote) bedside device.
 - 345 ○ The Cockpit will send an settings adjustment message to the (remote) bedside device.
 - The device will adjust the settings are requested; assuming the Cockpit is authorized to request the change.
 - The device will send a confirmation of the change to the Cockpit.

350

Post-conditions:

Reduced need for caregivers to put on and remove PPE and reduced likelihood of infection of caregivers.

355

X.4.2.4 Use Case #4: SDPi-A Silent PoC Alert Distribution Use Case - SPoC

Alerts from bedside devices are delegated to and signaled at a Central Station.

Additional information on this Use Case can be found in Volume 1 – Appendix C.5.

X.4.2.4.1 SDPi SPoC Use Case Description

360 This Use Case covers the implementation of a “Silent” Point of Care. SDC based devices in the patient room delegate their alerts to a multi-patient Central Station.

X.4.2.4.2 SDPi SPoC Process Flow

365 *<Diagram and describe the process flow(s) covered by this profile in order to satisfy the use cases. Demonstrate how the profile transactions are combined/sequenced. To provide context and demonstrate how the profile interacts with other profiles, feel free to include transactions and events that are “external” to this profile (using appropriate notation.)>*

The set of process flows will typically be exemplary, not exhaustive (i.e., it will address all the use cases, but will not show all possible combinations of actors, or all possible sequencing of transactions).

370 *If there are detailed behavioral rules that apply to a specific process flow or multiple process flows, an appendix may be added as needed.>*

<The roles at the top of the swimlane diagram should correspond to actor names, include the profile acronym:actor name if referencing an actor from a different profile.>

<Modify the following “Swimlane Diagram”.>

375




Figure X.4.2.4-1: Basic Process Flow in SDPi Profile – SpoC Use Case

380 *<If process flow “swimlane” diagrams require additional explanation to clarify conditional flows, or flow variations need to be described where alternate systems may be playing different actor roles, document those conditional flows here.>*

<Delete the material below if this is a workflow or transport profile. Delete the material above if this profile is a content module only profile.>

385 **Pre-conditions:**

- Bedside devices are SDC compliant
- At least one xDASmp (see Annex C) that is SDC compliant
- All devices in room delegate their alerts to one or more xDASmp (see Annex C).

390 **Main Flow:**

- Devices are all turned on.
 - Devices advertise their presence
- User sets the multi-patient Central Station (xDASmp) device into configuration mode.
 - Central Station device builds list of available devices of each room it is responsible for.
 - 395 ○ User configures which devices he/she wants to view on the Central Station.
 - Central Station advises viewed devices that it will be alerting on their behalf.
- User sets the Central Station device into operational mode.
 - Central Station connects to each configured device and collects the data from each device and displays it on the Central Station display.
 - 400 ○ Device goes into alarm
 - If a device viewed by the Central Station goes into alarm, the Central Station will signal the alarm.
 - The Central Station will acknowledge receipt and processing of the alarm back to the alarming device.
 - 405 ○ The alarming device may or may not signal the alert based on its configuration.
- User decides to Accept the alarm
 - If the user decides to accept the alarm, the Central Station will send an Accept message to the alarming device and stop signaling the alarm locally.

410

Post-conditions:

A patient room that has no audible alarms or alert sounds going off. This enables the patient to get more rest and less stressed.

415

X.4.2.5 Use Case #5: SDPi-A Silent ICU Alert Distribution - SICU

Alerts from bedside devices and Central Stations are delegated to and signaled at the caregiver's mobile device.

Additional information on this Use Case can be found in Volume 1 – Appendix C.6.

420 X.4.2.5.1 SDPi SICU Use Case Description

This Use Case covers the implementation of a “Silent” ICU. SDC based devices in the patient room delegate their alerts to a multi-patient Central Station which in turn delegates its alerts to a caregiver's mobile device. This results in a Silent ICU.

X.4.2.5.2 SDPi SICU Process Flow


425 *<Diagram and describe the process flow(s) covered by this profile in order to satisfy the use cases. Demonstrate how the profile transactions are combined/sequenced. To provide context and demonstrate how the profile interacts with other profiles, feel free to include transactions and events that are “external” to this profile (using appropriate notation.)>*

430 *The set of process flows will typically be exemplary, not exhaustive (i.e., it will address all the use cases, but will not show all possible combinations of actors, or all possible sequencing of transactions).*

If there are detailed behavioral rules that apply to a specific process flow or multiple process flows, an appendix may be added as needed.>

435 *<The roles at the top of the swimlane diagram should correspond to actor names, include the profile acronym:actor name if referencing an actor from a different profile.>*

<Modify the following “Swimlane Diagram”.>



<If process flow “swimlane” diagrams require additional explanation to clarify conditional flows, or flow variations need to be described where alternate systems may be playing different actor roles, document those conditional flows here.>

445 *<Delete the material below if this is a workflow or transport profile. Delete the material above if this profile is a content module only profile.>*

Pre-conditions:

- Bedside devices are SDC compliant
- An Alert Distribution System xDAScg (see Annex C) that is SDC compliant
- 450 • All Devices in Room delegate their alerts to one or more DASxx or CDASxx (see Annex C).

Main Flow:

- Devices are all turned on.
 - 455 ○ Devices advertise their presence
- User sets the caregiver Alert Distribution System (xDAScg) device into configuration mode.
 - The ADS device builds list of available devices of each room it is responsible for.
 - User configures which devices he/she wants to monitor on the ADS.
 - 460 ○ The ADS advises monitored devices that it will be alerting on their behalf.
 - The user also assigns specific rooms to specific caregivers.
- User sets the ADS device into operational mode.
 - The ADS connects to each configured device and monitors the alerts from each device.
- 465 • Device goes into alarm
 - If a device monitored by the ADS goes into alarm, the ADS will signal the alarm at the configured mobile device.
 - The ADS and mobile device will acknowledge receipt and processing of the alarm back to the alarming device.
 - 470 ○ The alarming device will not signal the alert. (It has delegated its alerts).

Post-conditions:

An ICU (or other Care Unit) that has no audible alarms or alert sounds going off. This enables the patients to get more rest and makes the overall Care Unit less stressful.

475

X.4.2.6 Use Case #6: SDPi-A/xC Remote Alert Management - RAM

Remote management of device alert setting using a personal device.

X.4.2.6.1 SDPi RAM Use Case Description

480 This Use Case covers the use of a personal device worn by an authorized caregiver to adjust alert related settings on a bedside device.

Additional information on this Use Case can be found in Volume 1 – Appendix C.7.

X.4.2.6.2 SDPi RAM Process Flow

485 *<Diagram and describe the process flow(s) covered by this profile in order to satisfy the use cases. Demonstrate how the profile transactions are combined/sequenced. To provide context and demonstrate how the profile interacts with other profiles, feel free to include transactions and events that are “external” to this profile (using appropriate notation.)>*

The set of process flows will typically be exemplary, not exhaustive (i.e., it will address all the use cases, but will not show all possible combinations of actors, or all possible sequencing of transactions).

490 *If there are detailed behavioral rules that apply to a specific process flow or multiple process flows, an appendix may be added as needed.>*

<The roles at the top of the swimlane diagram should correspond to actor names, include the profile acronym:actor name if referencing an actor from a different profile.>

<Modify the following “Swimlane Diagram”.>

495




Figure X.4.2.6-1: Basic Process Flow in SDPi Profile - RAM

500 *<If process flow “swimlane” diagrams require additional explanation to clarify conditional flows, or flow variations need to be described where alternate systems may be playing different actor roles, document those conditional flows here.>*

<Delete the material below if this is a workflow or transport profile. Delete the material above if this profile is a content module only profile.>

505 **Pre-conditions:**

- Bedside devices are SDC compliant
- A CDAScg (see Annex C) that is SDC compliant
- Devices in Room allow remote setting adjustment

510 **Main Flow:**

- User decides to adjust bedside setting from mobile device
 - The mobile device will send a settings adjustment message (via the ADS?) to the (remote) bedside device.
 - The device will adjust the settings are requested; assuming the mobile device is authorized to request the change.
 - The device will send a confirmation of the change to the mobile device.

515

Post-conditions:

Device settings are adjusted according to the caregiver’s input from a mobile device.

520

X.4.2.7 Use Case #7: SDPi-A Smart Alerting System - SAS

Integration of a Smart Alerting System into an ICU.

Additional information on this Use Case can be found in Volume 1 – Appendix C.8.

525 **X.4.2.7.1 SDPi SAS Use Case Description**

This Use Case covers the scenarios involved with the deployment of a separate Smart Alerting System in an ICU.

X.4.2.7.2 SDPi SAS Process Flow

530 *<Diagram and describe the process flow(s) covered by this profile in order to satisfy the use cases. Demonstrate how the profile transactions are combined/sequenced. To provide context*

and demonstrate how the profile interacts with other profiles, feel free to include transactions and events that are “external” to this profile (using appropriate notation.)

535 *The set of process flows will typically be exemplary, not exhaustive (i.e., it will address all the use cases, but will not show all possible combinations of actors, or all possible sequencing of transactions).*

If there are detailed behavioral rules that apply to a specific process flow or multiple process flows, an appendix may be added as needed.>

<The roles at the top of the swimlane diagram should correspond to actor names, include the profile acronym:actor name if referencing an actor from a different profile.>

540 *<Modify the following “Swimlane Diagram”.>*




Figure X.4.2.7-1: Basic Process Flow in <Profile Acronym> Profile

545 *<If process flow “swimlane” diagrams require additional explanation to clarify conditional flows, or flow variations need to be described where alternate systems may be playing different actor roles, document those conditional flows here.>*
<Delete the material below if this is a workflow or transport profile. Delete the material above if this profile is a content module only profile.>

550

Pre-conditions:

- Bedside devices are SDC compliant
- A single smart alert system that is SDC compliant
- The SAS is a component with no visual or audible alert outputs
- 555 • There is an xDASmp and/or xDAScg that is SDC compliant (see Annex C for explanation of these terms.)

Main Flow:

- Devices are all turned on.
 - 560 ○ Devices advertise their presence
- User sets the Smart Alert System (SAS) device into configuration mode.
 - The SAS builds a list of available devices in the rooms it is responsible for.
 - User configures which devices he/she wants to be monitored by the SAS.
 - The SAS advises viewed devices that it will be alerting on their behalf.
 - 565 ○ The SAS also is configured to connect to one or more Central Stations.
- User sets the SAS device into operational mode.
 - The SAS connects to each configured device and collects the data and alerts from each device.
 - If no alert is communicated from any of the configured devices, the SAS may still generate alert messages based on data trends or other factors it is monitoring. It communicates these Smart Alert to the Central Station(s).
 - 570
- Device goes into alarm
 - If an alert is communicated from any of the configured devices, the SAS will decide whether to forward the alert to the Central Station or create a Smart Alarm based on other data and alerts it has received.
 - 575
 - The SAS will acknowledge receipt and processing of the alarm back to the alarming device.

- The alarming device may or may not signal the alert locally based on its configuration.

580 *If the user decides to accept the alarm, the Cockpit will send an Accept message to the alarming device and stop signaling the alarm locally*

Post-conditions:

Smart Alerts are communicated properly with fallback modes in case of failures.

585

Appendices to Volume 1

<Add appendices to Volume 1 for this profile here. Examples of an appendix include HITSP mapping to IHE Use Cases or long use case definitions.>

<If there are no Volume 1 appendices, enter “Not applicable” and delete the Appendix A and Appendix B placeholder sections.>

590

<Volume 1 appendices are informational only. No “SHALL” language is allowed in a Volume 1 Appendix.>

Appendix C – Use Cases - Detail

595 This Appendix provides further details concerning the Use Cases that illustrate the use of the IHE SDPi actors and transactions.

C.1 Overview of the Concepts for DIS and DAS

600 While IEC 60601-1-8 is focused on alarm and alert functionality, it also provides some very useful system concepts such as the Distributed Information System (DIS) and Distributed Alarm System (DAS). We use these concepts in many of our use cases, so we have included the following is a quick guide to the functionality of DIS, CDIS, DAS and CDAS systems.

DIS – Distributed Information System

- DIS is a system for reporting alarm signals with no technical confirmation (of receipt).
 - 605 ○ Cannot rely on it for alarm signaling as a risk control
 - Optional support operator alarm management* response locally
 - Examples:
 - “Dashboard” - display which integrates the data from one patient. Dashboards do not support audible alerts or remote control.
 - 610 ▪ “View Station” - display which integrates the data from multiple patients. View Stations do not support audible alerts or remote control.

CDIS – Distributed Information System with Confirmation

- CDIS is a system for reporting alarm signals with no technical confirmation and operator confirmation (accept/reject). (Note it is not recognized in 60601-1-8)
 - 615 ○ Cannot rely on it for alarm signaling as a risk control
 - Optional support operator alarm management* response locally and remotely
 - Example – two-way pager (open loop)

xDIS – Connotes a system that can be either a DIS or CDIS

620

DAS – Distributed Alarm System

- DAS is a system for reporting alarm signals with technical confirmation (of receipt).
 - Can rely on it for alarm signaling as a risk control
 - Supports local alert confirmation
 - 625 • A communications failure or failure in any remote component of the DAS must initiate a technical alarm.
 - Examples:
 - “Cockpit”- display which integrates the data from one patient and support audible alerts. The DAS Cockpit only supports local operations while the

- 630 CDAS cockpit supports operator confirmation and optional remote control.
- “Central Station” - display which integrates the data from one patient and support audible alerts. The DAS Central Station only supports local operations while the CDAS Central Station supports operator confirmation and optional remote control.
- 635

CDAS - Distributed Information System with Confirmation

- CDAS is a system for reporting alarm signals with technical and operator confirmation (accept/reject) (of receipt).
- 640
- Can rely on it for alarm signaling as a risk control
 - Supports operator confirmation (accept/reject); It may redirect...
 - Optionally support local/remote alarm management (acknowledgement)
 - A communications failure or failure in any remote component of the DAS must initiate a technical alarm.
- 645
- Examples:
 - Cockpit and/or Central Station with remote confirmation and optional alarm management.
 - System that sends alarm to caregiver mobile device with accept / reject. Integrator may redirect

650

xDAS – Connotes a system that can be either a DAS or CDAS

In addition to the various types of xDIS and xDAS, we have also distinguished between various types of xDIS and xDAS as follows:

- 655
- xDxSsp – forms of DxS that are for single patients.
 - xDxSmp – forms of DxS that are for multiple patients
 - xDxScg – forms of DxS that communicated directly to the caregiver.

The following table provides examples and summarizes the various types of information/alarm systems.

660

Description	Type	Technical Delivery Confirmation ¹	Operator Delivery Confirmation ²	Optional Alarm Management	Examples
Reports alerts from a Single Patient (sp)	DISsp	No	No	Local	Single-Pt. information Dashboard
	CDISsp	No	Yes ³	Remote ³	Single-Pt. Remote View w/ accept/reject
	DASsp	Yes	No	Local	Single Pt. Cockpit w/audible alarms
	CDASsp	Yes	Yes	Remote	Single Pt. Cockpit w/ accept/reject
Reports alerts from Multiple Patients (mp)	DISmp	No	No	Local	Multiple-Pt. info. View Station
	CDISmp	No	Yes ³	Remote ³	Multiple-Pt. info. View Station w/ accept/reject
	DASmp	Yes	No	Local	Multiple Pt. Central Station w/ audible alarms
	CDASmp	Yes	Yes	Remote	Multiple Pt. Central Station w/ accept/reject
Reports and directs alerts to responsible Caregiver (cg)	DIScg	No	No	Local	Alert Distribution System (ADS) to caregiver pager, Mobile viewer
	CDIScg	No	Yes ³	Remote ³	ADS to caregiver pager, w/ accept/reject
	DAScg	Yes	No	Local	ADS to caregiver w/ audible/haptic alarms
	CDAScg	Yes	Yes	Remote	ADS to caregiver w/ accept/reject
¹ In each communication step the receiving device provides a technical response to the sending device that it received and is taking responsibility for the alert ² Operator can, at their choice, use the receiving device (communicator) UI to accept or reject responsibility for the alert ³ Not recommended since there is no confirmation that the Source has received the commands					

C.2 Feature: Use Case 1 - SDPi-R Standalone OR Dashboard - SORD

C.2.1 Narrative:

665 Dr. Presky is in the Operating Room performing an operation to remove cancerous tumors from
the patient's abdominal area. He can view previous radiology results, electrosurgical equipment
settings, patient readings such as HR, Blood Pressure, SpO2 and associated waveforms
integrated on his real-time 'Dashboard' display. The dashboard display can display visual
alarms but does not sound alerts or provide any remote-control capabilities. (This display can be
670 considered an xDISp.)

C.2.2 Background: Technical Pre-Conditions

Given all devices communicate using SDC

And at least one OR Dashboard display

675

C.2.3 Scenario: SORD 1.1 - OR Devices are Accessible to the Dashboard

Given dashboard detected at least one accessible OR device

When one or more OR Devices are operational

680 **Then** the Cockpit shall display parameter, waveform, setting, alarm, imaging, etc. information
from those devices (based on configuration)

C.2.4 Scenario: SORD 1.2 - OR Devices are inaccessible to the Dashboard

Given dashboard did not detect any accessible OR devices

Then the Dashboard shall display an error message

685

C.2.5 Scenario: SORD 1.3 - One or more OR Devices are inaccessible to the Dashboard

Given dashboard did not detect configured OR devices (based on configuration)

Then the Dashboard shall display an error message

690

C.3 Feature: Use Case 2 - SDPi-R/A/xC Standalone OR Cockpit - SORC

C.3.1 Narrative:

- 695 John Miller (13yrs, m) has chronic rhinosinusitis, which is an inflammatory condition in which the nose and his left maxillary sinus is swollen and the drainage of the mucus is prevented. John's chronic rhinosinusitis doesn't respond to medication anymore. After consulting with his physician, he and his parents decide to resolve the issue with Functional Endoscopic Sinus Surgery (FESS). The FESS will be done as a day surgery, so that John can get home in the evening.
- 700 Before the day of the surgery, a CT scan is taken that is used to guide the surgeon during the surgery.
- In order for the surgery to start, John is put under general anesthesia and monitored with a patient monitor by a pediatric anesthesiologist, especially his mean arterial blood pressure which has been lowered in order reduce capillary bleeding to provide optimal visibility of the surgical field.
- 705 During the intervention, the Surgeon has a constant view, using his Surgical Cockpit of the patient's vitals (including MABP), CT imaging results, real-time endoscope camera output and has access to the control functions to execute the intervention. The anesthesiologist can also view relevant patient real-time information such as ECG, blood pressure, anesthesia agent, depth of anesthesia, allergies, etc. using the Anesthesia Cockpit where he/she can manage alarms and
- 710 control device settings as needed.
- {During the procedure one of the surgical devices has a technical issue. It generates a technical alert which notifies the responsible biomedical technician. He/she decides to replace the device and connects it to the network where it is automatically discovered and configured allowing the intervention to continue.}
- 715 There are no additional technical or clinical problems, the surgery is a success and John can go home with his parents.

C.3.2 Background: Technical Pre-Condition

Given all devices communicate using SDC

- 720 **And** at least one Anesthesia Cockpit and/or Surgical Cockpit

C.3.3 Scenario: SORC 2.1 - OR Devices are accessible to the Cockpit

Given Cockpit detected at least one accessible OR device

When the one or more OR Devices are operational

- 725 **And** audio alarms are enabled on the Cockpit

And remote control is enabled on the Cockpit

Then the Cockpit shall display parameter, waveform, setting, alarm, etc. information from those devices (based on configuration)

730 **And** the audio alarm shall be disabled on the source devices associated with the Cockpit (based on configuration)

And the user shall be able to control device settings at the source device or at the Cockpit

C.3.4 Scenario: SORC 2.2 - OR Devices are inaccessible to the Cockpit

Given Cockpit did not detect any accessible OR devices

735 **When** any OR Device detects an alert condition

Then the Cockpit will display an error message

And the audio alarm shall be signaled at that source device in the OR

C.3.5 Scenario: SORC 2.3 - Some OR Devices are inaccessible to the Cockpit

740 **Given** Cockpit did not detect some configured OR devices (based on configuration)

When the OR Device detects an alert condition

Then the Cockpit will display an error message

And the audio alarm shall be signaled at that source device in the OR for the devices not detected by the cockpit

745

C.4 Feature: Use Case 3 - SDPI-R/xC ICU Isolation PoC Use Case - IIPoC

Narrative:


750 **Background: Technical Pre-Condition**

Coming soon...

755 **C.5 - Feature: Use Case 4 - SDPI-A - Silent PoC Alert distribution - SPoC**

C.5.1 Narrative

Linda is an ICU nurse responsible for 4 patients. While she is updating documentation at the nursing station Jack's (one of her patients) condition deteriorates and his ventilator goes into an alarm state (#1). The ventilator alarm sounds are quite loud and jarring which usually disturbs the patient in the room and nearby rooms. In this case the alert only generates visual alerts and does not generate an audible tone in the patient room (#2). It does generate an alarm tone at the central station (xDASmp) (#3) and her mobile device (xDAScg) (#4). As a result, Linda must acknowledge or otherwise handle the alert at the bedside (#5).



C.5.2 Background: Technical Pre-Condition

Given all devices communicate using SDC

And at least one xDASmp

And all devices in room delegate their alerts to one or more xDASmp.

775 **C.5.3 Scenario: SPoC 4.1 - Medical device detects an alert situation and at least one distributed alarm system (xDASxx) is accessible**

Given alert event was detected by a medical device attached to the patient

When at least one remote alert system is accessible

Then the alert shall be shown on all accessible remote alerting devices

And the audio alarm shall be enabled on all accessible remote alerting devices

780 **And** active device alert events shall be shown on the medical devices locally

And the audio alarm shall be disabled on all medical devices in the patient room

785 **C.5.4 Scenario: SPoC 4.2 - Medical device detects an alert situation and all distributed alarm systems (xDASxx) are inaccessible or become inaccessible**

Given alert event was detected by a medical device attached to the patient

When distributed alarm systems (DAS/DIS) are inaccessible or become inaccessible

Then active device alert events shall be shown on the medical devices locally

And the audio alarm shall be enabled on all medical devices in the patient room

790

C.5.5 Scenario: SPoC 4.3 - Alert situation has been resolved and at least one distributed alarm system (xDASxx) is accessible

Given medical device detected that the alert situation has been resolved

When at least one distributed alarm system (DAS/DIS) is accessible

795 **Then** the alert shall be shown as inactive/ended at the medical device locally

And the audio alarm shall be disabled on the medical device in the patient room

And the alert shall be shown as inactive/ended on all accessible remote alerting devices


And the audio alarm shall be disabled on all accessible remote alerting device for this alert event

800 **C.5.6 Scenario: SPoC 4.4 - Medical device detects an alert situation, initially xDASxx is accessible but fails**

C.6 Feature: Use Case 5 - SDPi-A: Silent ICU Alert Distribution - SICU

805 **C.6.1 Narrative:**

810 Jackie is an ICU nurse responsible for 4 patients. While she is updating documentation at the nursing station Sven's (one of her patients) condition deteriorates and his ventilator goes into an alarm state (#1). The ventilator alarm sounds are quite loud and jarring which usually disturbs the patient in the room and nearby rooms. Jackie's ICU uses personal mobile devices to alert the nurses about patient alarms so in this case the alert only generates visual alerts in the patient room (#2) and central station (#3). It does generate an audible alert on her mobile device (xDAScsg) (#4).



815

C.6.2 Background: Technical Pre-Condition

Given all devices communicate using SDC

And at least one xDAScsg

And all devices in room delegate their alerts to one or more xDASxx

820

C.6.3 Scenario: SICU 5.1 - Medical device detects an alert situation and the distributed alarm system CDAScsg is accessible

Given alert event was detected by a medical device attached to the patient

When at least one remote alert system is accessible

825 **Then** the alert shall be shown on all accessible remote alerting devices

And the audio alarm shall be enabled on the caregiver's accessible CDAScsg remote alerting devices

And the audio alarm shall be disabled on all accessible non-CDAScsg remote alerting devices

And active device alert events shall be shown on the medical devices locally

830 **And** the audio alarm shall be disabled on all medical devices in the patient room

C.6.4 Scenario: SICU 5.2 - Medical device detects an alert situation and the distributed alarm system CDAScsg is inaccessible

Given alert event was detected by a medical device attached to the patient

835 **When** at least one remote alert system is accessible

Then the alert shall be shown on all accessible remote alerting devices

And the audio alarm shall be enabled on all accessible remote alerting devices

And active device alert events shall be shown on the medical devices locally

And the audio alarm shall be disabled on all medical devices in the patient room

840

C.6.5 Scenario: SICU 5.3 - Medical device detects an alert situation and all distributed alarm systems (xDASxx) are or become inaccessible

Given alert event was detected by a medical device attached to the patient

When distributed alarm systems (DAS/DIS) are inaccessible or become inaccessible

845 **Then** active device alert events shall be shown on the medical devices locally
And the audio alarm shall be enabled on all medical devices in the patient room

C.6.6 Scenario: SICU 5.4 - Alert situation has been resolved and at least one distributed alarm system (xDASxx) is accessible

850 **Given** medical device detected that the alert situation has been resolved
When at least one distributed alarm system (DAS/DIS) is accessible
Then the alert shall be shown as inactive/ended at the medical device locally
And the audio alarm shall be disabled on the medical device in the patient room
And the alert shall be shown as inactive/ended on all accessible remote alerting devices
855 **And** the audio alarm shall be disabled on all accessible remote alerting device for this alert event


C.6.7 Scenario: SICU 5.5 - Medical device detects an alert situation, initially DAS is accessible but fails

860

C.7 Feature: Use Case 6 - SDPi-A/xC: Remote Alert Management - RAM

C.7.1 Narrative

865 Boris is an ICU nurse responsible for 4 patients. His ICU has a central station but also uses personal devices for alert notification and management. He needs to adjust the upper heart rate limit for Caroline, one of his patients. Even though Boris is near the central station, he decides to use his personal device to adjust the limit (#1).



870

C.7.2 Background: Technical Pre-Condition

Given all devices communicate using SDC

875 **And** at least one CDAScg

And all devices in room enable remote control

C.7.3 Scenario: RAM 6.1 - Caregiver adjusts alarm limit at their Mobile Device

Given alert event was detected by a medical device attached to the patient

880 **And** remote alerting device is part of the CDAScg

When caregiver confirms the alert at a remote alerting device

And CDAScg is accessible

Then the alert shall be shown as acknowledged at the medical device

And the audio alarm shall be disabled on the medical device

885 **And** the alert shall be shown as acknowledged on all accessible remote alerting devices

C.7.4 Scenario: RAM 6.2 - Caregiver attempts to adjust alarm limit at their CDAScg mobile device but the respective CDAScg is not accessible

Given alert event was detected by a medical device attached to the patient

890 **And** remote alerting device is part of the DIS

When caregiver confirms the alert at a remote alerting device

But DAS is inaccessible

Then ???

895

C.8 Feature: Use Case 7 - SDPi-A Smart Alerting System - SAS


C.8.1 Narrative:

900 Ben is an ICU nurse. Suddenly, he gets a “Check Ventilation Hose!” alert for one of the patients he is responsible for. In addition, the alert suggests possible root causes to the caregiver (i.e. obstruction (sputum/kinked hose)). This leads to a quicker, more adequate intervention.

905 The alert was generated by a smart alerting system that collects all the data from the point-of-care devices such as vital signs, alerts, settings, waves, etc., and combines them to create more actionable information for the care giver to guide care, intervention and treatment. In the example above, an algorithm combines a “Low SPO2” alarm from the patient monitor and a “Peak Pressure” alarm and “Minute Volume low” alarm from the ventilator into one alarm superseding the individual alarms.

910 The original alerts generated by the patient monitor and the ventilator are shown at the devices but the audio alarm is enabled or disabled on both devices dependent on other rules such as configuration or presence of caregiver in patient room.

915 Note that the smart alerting system is seen as a separate entity independent of an xDxSxx in this Feature. However, a combination of a xDxSxx with a smart alerting system is a typical configuration.



920

C.8.2 Background: Technical Pre-Conditions

Given bedside devices are SDC compliant

And a single smart alert system (SAS) that is SDC compliant

And the SAS is a component with no visual or audible alert output

925 **And** an xDASmp and/or xDAScg that is SDC compliant

C.8.3 Scenario Outline: SAS 7.1 - Local device generates alerts, Smart Alerting System is accessible, and xDASmp and/or xDAScg is accessible

Given local device audio alarm state was set to <state>

930 **When** there is an alert event on one or more medical devices in the patient room

And smart alerting system is accessible

And xDASmp and/or xDAScg is accessible

Then the alerts on the medical devices in the patient room shall be delegated to the SAS

And active device alert events shall be shown on the medical devices locally

935 **And** the audio alarm shall be <action> on all medical devices in the patient room

And smart alerts from the SAS shall be delegated to the xDASmp and/or xDAScg

Examples:

	state	action	
940	disabled	disabled	
	enabled	enabled	

C.8.3 Scenario Outline: SAS 7.2 - Local device generates alerts, Smart Alerting System is inaccessible, and xDASmp and/or xDAScg is accessible

945 **Given** local device audio alarm state was set to <state>

When there is an alert event on one or more medical devices in the patient room

And smart alerting system is inaccessible

And xDASmp and/or xDAScg is accessible

950 **Then** the alerts on the medical devices in the patient room shall be delegated to the xDASmp and/or xDAScg

And active device alert events shall be shown on the medical devices locally

And the audio alarm shall be <action> on all medical devices in the patient room

Examples:

955	state	action	
	disabled	disabled	
	enabled	enabled	

C.8.4 Scenario: SAS 7.3 - Local device generates alerts, and xDASmp and/or xDAScg is inaccessible

960 **When** there is an alert event on one or more medical devices in the patient room

And xDASmp and/or xDAScg is inaccessible

Then active device alert events shall be shown on the medical devices locally

And the audio alarm shall be enabled on all medical devices in the patient room

965

C.8.5 Scenario: SAS 7.4 -Smart Alerting generates an alert and an xDASxx is accessible

When SAS is accessible

And there is an alert event detected by the SAS (e.g. derived from vital signs data)

970 **And** an xDASxx is accessible

Then the alerts on the SAS shall be delegated to the xDASxx

And the xDASxx shall signal the SAS audio and visual alerts

And a CDASxx can confirm the alert back to the SAS

975 **C.8.6 Scenario: SAS 7.5 - Smart Alerting generates an alert and an xDASxx is inaccessible**

When SAS is accessible

And there is an alert event detected by the SAS

And an xDASxx is inaccessible

980 **Then** any device alert signal delegation shall be disabled

And any active alerts (audio and visual) shall be signaled on the originating device

C.8.6 Scenario: SAS 7.6 - Smart Alerting System is inaccessible from devices and an xDASxx is accessible

985 **When** SAS is inaccessible

And an alert event was detected by a bedside device

And an xDASxx is accessible

Then the alerts on the device shall be delegated to the xDASxx

And the xDASxx shall signal the device audio and visual alerts

990 **And** a CDASxx can confirm the alert back to the device

C.8.7 Scenario: SAS 7.7 - Smart Alerting System is inaccessible from devices and an xDASxx is inaccessible

When SAS is inaccessible

995 **And** an alert event was detected by a bedside device

And an xDASxx is inaccessible

Then any device alert signal delegation shall be disabled

And any active alerts (audio and visual) shall be signaled on the originating device

1000 **C.8.8 SAS Scenario 7.8: Devices, CDASxx and Smart Alerting System are accessible, and SAS generates an alert**

When a CDASxx is accessible

And an alert event was detected by the SAS

And the operator confirms the alert at the CDASxx

1005 **Then** the confirmation will be sent to the SAS

And the SAS will send the confirmation to the appropriate Devices